



Maximum Load: One Person

SWL: 260kg

For both vertical and horizontal lifting, based on a Factor of Safety of 10:1

Fitting and adjustment

Undo all straps.

Lay casualty on stretcher pass leg strap between the legs and attach shoulder straps (black with orange fleck) .

Fasten all black straps - tighten snugly.

Fasten head strap around head if necessary.

When lift commences and the casualty's weight shifts adjust all straps if necessary for a secure fit.

In use

Carry the stretcher using any of the handles along the sides. DO NOT carry it by the red PVC loop at the head end as this may flex the casualty's neck.

The stretcher does not provide full spinal immobilisation. Use separate splinting for spinal or limb injuries as necessary.

Vertical Lift

Attach all ropes to the red PVC loop at the head end, this loads the stretcher and integral harness simultaneously.

Horizontal Lift

Attach the lifting slings provided to the side handles with a 'choke' hitch. The lifting slings are labelled "HEAD", "MIDDLE" and "FEET", according to their position on the stretcher. Connect all the slings to the anchor plate. Attach all ropes to the anchor plate or if required the backup rope may be attached to the red PVC loop at the head end, this will connect to the casualty via the integral harness.



